



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

ZAP™ SLIK 200 uL Low Retention Aerosol Filter Pipet Tips for Rainin® LTS Style Pipettors, in 96 Racks, Sterile

Fluid Measure	μL	Freezable	NA
Fluid Volume	200	Freeze Temp	
Height (mm)	51.3	Autoclavable	121°C 15min
Width (mm)		Sterile	Yes
Length (mm)		Sterile Processing	Certified as sterile with a sterility assurance level (SAL) of 10 ⁻⁶
Pack Qty	960	Centrifugation	NA
Case Qty	9600	Centrifugation Speed (RCF)	
Package Type	96 Rack	Shelf Life (Months)	36
Material	Proprietary Polypropylene	Certification	
Cap Material		Regulatory	

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.

Additional Testing: Protease, Human DNA, ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)



This Technical Data Sheet is provided as a service and for documentation only. Labcon North America retains all rights to the information, procedures, and results of the material contained hereon as proprietary information. Any copying and distribution of this technical data without prior written authorization is prohibited by law. The use of this information and documentation does not discharge the user from their responsibility to test and use this product appropriately. This item is for use solely by laboratory professionals. This document is automatically generated and is © **Copyright 10-May-2021** . All Rights Reserved. Technical Data Sheet Version 1/2020